Modern ergometer with multifunctional applications







Highlights

High standards

Lode is a socially and environmentally responsible company. All Lode products are RoHS/WEE compliant and Lode is ISO 9001:2003, ISO 13485:2008 and FDA 510K certified. All medical products comply to MDD 93/42/EEC, incl. IEC 60601-1.

Easy to operate

For Lode products this means:

- easy to connect
- easy to move around
- easy user interface

Reliable and reproducible stress tests

The experience of professionals who calibrate many ergometers shows that the Lode ergometers are the most reliable across the complete workload and rpm range and still within specifications even after many years of intensive use.



Multifunctional

The ergometer can be used in various ergometry settings, enabling a multifunctional deployment.

Additional features with PCU

Besides the possibility to program 24 protocols easily, this control unit offers the following features:

- better monitoring because of the additional and larger display

- a perfect combination with BPM
- possibility to measure SpO2



Modern ergometer with multifunctional applications



The Angio cpet is an ergometer that can be used for both arm and supine ergometry. Its compact design makes it universally applicable for ergometric studies in those sectors in which standard ergometry cannot be used. The Angio operates independent of pedaling speed in the range of 7 - 1000 watt. The Angio cpet is standard supplied with a communication module and can therefor be easily controlled by all known stress ECG and pulmonary devices in the world. The workload, rpm and time can be readout from the 3,5" colour display. The ergometer is supplied with standard pedals.

Features



Accurate over a long period of time

The Lode ergometers are supplied with an electro-magnetic braking mechanism of Lanooy (eddy current). The biggest advantage of this braking system compared to a friction braking system is the absolute

to a friction braking system is the absolute accuracy and the accuracy over time. Moreover, friction braking systems have more wearing parts.



Extreme low start-up load

The extreme low start-up load of 7 watts and the adjustability in small steps of 1 watt make this ergometer perfectly suitable for many different applications. The standard control unit shows multiple ergometry parameters and you can determine your specific default setting and start-up menu.



Small adjustment steps

The workload of the Lode ergometers is adjustable in steps of only 1 watt. Depending on your wishes, the test operator or the test subject can adjust the workload. The steps of 1 watt are possible in the manual mode as well as within protocols.



LEM compatible

This product can be used with Lode Ergometry Manager (LEM) software to manage data and to apply specific protocols when a Communication card or Network card is present



Versatile

controls

LCRM compatible

This product can be used with Lode Cardiac Rehabilitation Manager software (LCRM)

Additional features with PCU

Besides the possibility to program 24 protocols easily, this control unit offers the following features:

- better monitoring because of the additional and larger display
- a perfect combination with BPM
- possibility to measure SpO2

Versatile Interfacing



Various interface protocols guarantee perfect communication with all commonly known stress ECG and spirometry equipment



Modern ergometer with multifunctional applications



Angio cpet - single set can a.o be extended with the following options:

Control Unit with touch screen 7" for ergometer Multifunctionality	Programmable Control Unit with 7" Touchscreen Programmable	Blood Pressure with ECG trigger for bicycle ergometer with ECG trigger	Heart rate for bicycle ergometers Heart rate in beats per minute	SpO2 for control unit with touch panel (bicycle) Saturation and heart rate
Partnumber: 945834	Partnumber: 945835	Partnumber: 945828	Partnumber: 945821	Partnumber: 945823
Pedal shoes (pair)	Pedal shoes pediatric (pair)	Pedal shoes extra large (pair)	Adjustable cranks	Handgrips (pair)
Extra stability during cycling	Pedal shoes for childen	For large feet sizes	Optimal force application	Versatile ergometry
			and the second	
Partnumber: 917803	Partnumber: 917833	Partnumber: 917834	Partnumber: 928804	Partnumber: 917812
Handgrip fixation set	USB to Serial converter	Network Module		
Handgrip fixation set	Easy connection	Easy networking with LEM and LCRM		
Partnumber: 917822	Partnumber: 226012	Partnumber: 945851		



Modern ergometer with multifunctional applications



Specifications

Workload				
Minimum load	7 W			
Maximum peak load	1000 W			
Minimum load increments	1 W			
Maximum continuous load	750 W			
Hyperbolic workload control	\checkmark			
Linear workload control	\checkmark			
Fixed torque workload control	\checkmark			
Maximum rpm independent constant load	150 rpm			
Minimum rpm independent constant load	30 rpm			
Optional heart rate controlled workload	\checkmark			
Electromagnetic "eddy current" braking system	\checkmark			
Accuracy				
Workload accuracy from 7 to 100 W	3 W			
Workload accuracy from 100 to 500 W	3 %			
Workload accuracy from 500 to 1000 W	5 %			
User Interface				
Readout Distance	\checkmark			
Readout RPM	\checkmark			
Readout Heartrate	\checkmark			
Readout target HR	\checkmark			
Readout Energy	\checkmark			
ReadoutTorque	\checkmark			
Readout Time	\checkmark			
Readout Power	\checkmark			
Set Display	\checkmark			
Set Resistance	\checkmark			
Set P-Slope	\checkmark			
Set Mode	\checkmark			
Manual operation mode	\checkmark			
Preset protocol operation mode	\checkmark			
External control unit	\checkmark			
Selfdesigned protocol operation mode	\checkmark			
Dimensions				
Productweight	26 kg	57.3 lbs		

Power requirements	
VAC	100-240 V
Phases	1
Frequency	50/60 Hz
Power consumption	160 W
Power cord IEC 60320 C13 with CEE 7/7 plug	\checkmark
Power cord NEMA	×
Standards & Safety	
IEC 60601-1:2012	\checkmark
ISO 13485:2003 compliant	\checkmark
ISO 9001:2008 compliant	\checkmark
Certification	
CTüVus according to NRTL - pending	\checkmark
CB according to IECEE CB - pending	\checkmark
CE class Im according to MDD93/42/EEC - pending	\checkmark

Order info

Partnumber:

r: 967900

*Specifications are subject to change without notice.



Distributed by 29 Costa Rica 7 Costa Rica 7 Tel: 6

Lode B.V. Zernikepark 16 9747 AN Groningen The Netherlands Tel: +31 50 5712811 Fax: +31 50 5716746 E-mail: ask@lode.nl Internet: www.lode.nl